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## Dafilon®

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<b>Instructions for use</b>	<b>GB</b>
<b>Dafilon®</b>	
<b>Description</b>	
Dafilon® is a sterile, monofilament, non-absorbable surgical suture material produced from polyamide. Dafilon® is available in blue (polyamide 6/66), coloured with Copper Phthalocyanine (C.I. 74 160) or in black (polyamide 6/6), coloured with Hematein Logwood Black C.I. 75290. To make it more readily recognizable, but is also available undyed in the natural beige colour. Dafilon® fulfils all the requirements of the European Pharm. and the United States Pharm. current edition for sterile, nonabsorbable sutures.	
<b>Indications</b>	

Dafilon® is basically indicated for dermal sutures, Microsurgery and Ophthalmology for ophthalmic, microvascular and neural interventions.

**Mode of action**

Suture materials are used primarily for adaptation of the wound edges to render possible an undisturbed wound healing. When Dafilon® suture materials are employed, a minimal inflammatory tissue reaction will occur followed by a gradual encapsulation of the suture material by fibrous connective tissue. Like all polyamide materials, over a long period of time Dafilon® in vivo is subject to a gradual loss of the total suture strength through the action of body enzymes. When Dafilon is used for skin closure, and being a non absorbable material, it is strongly recommended to remove the suture once the healing period is achieved, and no later than 30 days.

**Contra-indications**

On account of the gradual loss of strength, which can occur in vivo over a longer period, the use of Dafilon® is contraindicated when a permanent retention of suture strength is of importance. Dafilon® microsutures having a limited tensile strength retention are not recommended for conventional woundclosure in intracutaneous/subcutaneous, cardiovascular and gastro-intestinal indications as well as for the fixation of vascular prostheses.

**Warning note**

Dafilon® must not be resterilized. Opened unused or damaged packs shall be discarded.

Do not reuse: Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product.

No special storage conditions are required.

**Note/precautionary measures**

When working with suture materials great care should be taken to ensure that the use of surgical instruments, such as tweezers and needle holders, does not cause the material to be damaged by being pinched or kinked. The user should be familiar with surgical suturing techniques, before employing Dafilon® suture materials.

**Application**

To be used in accordance with surgical requirements.

**Side effects**

Unfavourable reactions are not known to occur when the material is employed according to the instructions.

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