General Information

3M™ Filtek™ Supreme XTE Flowable Restorative, is a low viscosity, visible-light activated, radiopaque, flowable nanocomposite. The restorative is packaged in capsules and syringes. It is available in a variety of tooth-colored shades. The shade offered with Filtek Supreme XTE flowable restorative were designed to coordinate with shades offered with Filtek™ Supreme XTE Universal Restorative, Filtek Supreme XTE flowable restorative contains bisGMA, TEGDMA and Proceryl resins. The fillers are a combination of yttrium trifluoride filler with a range of particle sizes from 0.1 to 5.0 microns, a non-agglomerated/non-aggregated surface-modified 20nm silica filler, a surface-modified aggregated zirconia/silica cluster filler (comprised of 20 nm silica and 4 to 11 nm zirconia particles). The aggregate has an average cluster particle size of 0.6 to 10 microns. The inorganic filler loading is approximately 65% by weight (46% by volume).

Indications

- Class III and V restorations
- Restoration of minimally invasive cavity preparations (including small, non stress-bearing occlusal restorations)
- Base/liner under direct restorations
- Repair of small defects in esthetic indirect restorations

Pit and fissure sealant**Undercut blockout****Repair of resin and acrylic temporary materials****Precautionary Information for Patients**

This product contains substances that may cause an allergic reaction by skin contact in certain individuals. Avoid use of this product in patients with known acrylate allergies. If prolonged contact with oral soft tissue occurs, flush with large amounts of water. If allergic reaction occurs, seek medical attention as needed. Remove the product if necessary and discontinue future use of the product.

Precautionary Information for Dental Personnel

This product contains substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured product. If skin contact occurs, wash skin with soap and water. Use of protective gloves and a no-touch technique is recommended. Avoid skin contact with uncured product. If product contacts skin, remove and discard glove, wash hands immediately with soap and water and then re-glove. If allergic reaction occurs, seek medical attention as needed.

3M MSDS information can be obtained from 3M.com/dental or contact your local subsidiary.

Instructions for Use**Preparation**

- Prophy: Teeth should be cleaned with pumice and water to remove surface stains.
- Shade Selection: Prior to isolation of tooth, select the appropriate shade(s) of Filtek Supreme XTE flowable restorative using a standard VITAPAN® classical shade guide. Shade selection accuracy may be enhanced via a mock-up; place the chosen shade of the tooth on the shade guide and compare the shade of the tooth with the shade guide.
- Isolation: A rubber dam is the preferred method of isolation. Cotton rolls and an evacuator can also be used.

Directions

- Protection: If a pulp exposure has occurred and the situation warrants a direct pulp capping procedure, use a minimum amount of calcium hydroxide on the exposure followed by an application of 3M™ Vitrebond® or Vitrebond® Plus Light Cure Glass Ionomer, Vitrebond or Vitrebond Plus liner/bases may also be used to limit area of deep cavity excavation.

2. Adhesive System: Follow manufacturer's instructions regarding placement.

Note: Silane treatment is recommended for repair of ceramic restorations followed by the adhesive application.

3. Dispensing:

- Syringe delivery: Filtek Supreme XTE flowable restorative can be delivered directly from the dispensing tip.

3.1. Protective eyewear for patients and staff is recommended when using the dispensing tip.

3.2. Prepare the dispensing tip: Remove cap and save. Twist the prebondable dispensing tip securely onto the syringe. Holding the tip away from the patient and any dental staff, express a small amount of Filtek Supreme XTE flowable restorative to assure that the delivery system is not plugged.

3.3. If plugged, remove the dispensing tip and express a small amount of Filtek Supreme XTE flowable restorative directly from the dispensing tip. If present, from the syringe opening. Replace dispensing tip and again express composite. Filtek Supreme XTE flowable restorative may be extruded onto a dispensing pad and applied with a brush or other appropriate instrument.

3.2. Single-Dose Capsule: Insert capsule into 3M™ Restorative Dispenser. Refer to separate restorative dispenser instructions for full instructions and precautions. Extrude restorative directly into cavity.

4. Placement:

4.1. Place and light cure restorative in increments as indicated in Section 5.

4.2. Avoid intense light in the working field.

4.3. To prevent oozing of paste from the syringe tip after placement, pull back on syringe plunger.

5. Curing: This product is intended to be cured by exposure to a halogen or LED light with a minimum intensity of 400 mW/cm² in the 400-500nm range. Cure each increment by exposing its entire surface to a high intensity visible light source, such as a 3M curing light. Hold the light guide tip as close to the restorative as possible during light exposure.

6. Complete the restoration:

6.1. Base/Liner application: Place a composite restorative material, such as Filtek Supreme XTE Universal Restorative directly over the cured Filtek Supreme XTE flowable restorative. Follow the manufacturer's instructions regarding placement, curing, finishing, occlusal adjustment, and polishing.

6.2. Direct Restorative Application:

6.2.1. Contour restoration surfaces with fine finishing diamonds, burs or stones.

6.2.2. Check occlusion with an articulating paper. Examine centric and lateral excursion contacts. Carefully adjust occlusion by removing material with a fine polishing diamond or stone.

6.2.3. Polish with 3M™ Sof-Lex™ Finishing and Polishing System or with white stones, rubber points or polished paste where discs are not suitable.

6.3. Pit and Fissure Sealants: Gently remove the inhibited layer remaining after light curing with slurry of pumice or polishing paste.

Storage and Use

1. This product is designed to be used at room temperature. If stored in cooler allow product to reach room temperature prior to use.

2. Shelf life at room temperature is 24 months. Ambient temperatures routinely higher than 27°C/80°F may reduce shelf life. See outer package for expiration date.

- Do not expose restorative materials to elevated temperatures, or to intense light.
- Do not store materials in proximity to eugenol containing products.
- Discard used syringe tip and replace with syringe storage cap.

Disinfect the products using an intermediate level disinfection process (liquid contact) as recommended by the Centers for Disease Control and endorsed by the American Dental Association, Guidelines for Infection Control in Dental Health-Care Settings, 2003. MMWR, December 19, 2003:52(RR-17), Centers for Disease Control and Prevention.

Disposal: See the Material Safety Data Sheet (available at 3M.com/dental or through your local subsidiary) for disposal information.

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Caution: U.S. Federal Law restricts this device to sale or use on the order of a dental professional.

Warranty

3M warrants this product will be free from defects in material and manufacture. 3M MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M's sole obligation shall be repair or replacement of the 3M product.

Limitation of Liability

Except where prohibited by law, 3M will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, notwithstanding the theory asserted, including warranty, contract, negligence or strict liability.

Precautionary Information for Patients

This product contains substances that may cause an allergic reaction by skin contact in certain individuals. Avoid use of this product in patients with known acrylate allergies. If prolonged contact with oral soft tissue occurs, flush with large amounts of water. If allergic reaction occurs, seek medical attention as needed. Remove the product if necessary and discontinue future use of the product.

Precautionary Information for Dental Personnel

This product contains substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured product. If skin contact occurs, wash skin with soap and water. Use of protective gloves and a no-touch technique is recommended. Avoid skin contact with uncured product. If product contacts skin, remove and discard glove, wash hands immediately with soap and water and then re-glove. If allergic reaction occurs, seek medical attention as needed.

3M MSDS information can be obtained from 3M.com/dental or contact your local subsidiary.

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3. Isolation: A rubber dam is the preferred method of isolation.

Cotton rolls and an evacuator can also be used.

4. Directions:

- Protection: If a pulp exposure has occurred and the situation warrants a direct pulp capping procedure, use a minimum amount of calcium hydroxide on the exposure followed by an application of 3M™ Vitrebond® or Vitrebond® Plus Light Cure Glass Ionomer, Vitrebond or Vitrebond Plus liner/bases may also be used to limit area of deep cavity excavation.

2. Adhesive System: Follow manufacturer's instructions regarding placement.

Note: Silane treatment is recommended for repair of ceramic restorations followed by the adhesive application.

3. Dispensing:

- Syringe delivery: Filtek Supreme XTE flowable restorative can be delivered directly from the dispensing tip.

3.1. Protective eyewear for patients and staff is recommended when using the dispensing tip.

3.2. Prepare the dispensing tip: Remove cap and save. Twist the prebondable dispensing tip securely onto the syringe. Holding the tip away from the patient and any dental staff, express a small amount of Filtek Supreme XTE flowable restorative to assure that the delivery system is not plugged.

3.3. If plugged, remove the dispensing tip and express a small amount of Filtek Supreme XTE flowable restorative directly from the dispensing tip. If present, from the syringe opening. Replace dispensing tip and again express composite. Filtek Supreme XTE flowable restorative may be extruded onto a dispensing pad and applied with a brush or other appropriate instrument.

3.2. Single-Dose Capsule: Insert capsule into 3M™ Restorative Dispenser. Refer to separate restorative dispenser instructions for full instructions and precautions. Extrude restorative directly into cavity.**4. Placement:**

4.1. Place and light cure restorative in increments as indicated in Section 5.

4.2. Avoid intense light in the working field.

4.3. To prevent oozing of paste from the syringe tip after placement, pull back on syringe plunger.

5. Curing: This product is intended to be cured by exposure to a halogen or LED light with a minimum intensity of 400 mW/cm² in the 400-500nm range. Cure each increment by exposing its entire surface to a high intensity visible light source, such as a 3M curing light. Hold the light guide tip as close to the restorative as possible during light exposure.

6. Complete the restoration:

6.1. Base/Liner application: Place a composite restorative material, such as Filtek Supreme XTE Universal Restorative directly over the cured Filtek Supreme XTE flowable restorative. Follow the manufacturer's instructions regarding placement, curing, finishing, occlusal adjustment, and polishing.

6.2. Direct Restorative Application:

6.2.1. Contour restoration surfaces with fine finishing diamonds, burs or stones.

6.2.2. Check occlusion with an articulating paper. Examine centric and lateral excursion contacts. Carefully adjust occlusion by removing material with a fine polishing diamond or stone.

6.2.3. Polish with 3M™ Sof-Lex™ Finishing and Polishing System or with white stones, rubber points or polished paste where discs are not suitable.

Storage and Use

1. This product is designed to be used at room temperature. If stored in cooler allow product to reach room temperature prior to use.

2. Shelf life at room temperature is 24 months. Ambient temperatures routinely higher than 27°C/80°F may reduce shelf life. See outer package for expiration date.

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