3M ESPE

Impregum™ Penta® Soft Quick

Ⓔ Polyether Impression Material

Ⓓ Polyether Abformmasse

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Brugsanvisning

Brugsanvisning

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Dental Products
Carl-Zeiss-Str. 1
41453 Neuss – Germany

3M, ESPE, Impregum, Penta, PentaMatric and Pentamix are trademarks of 3M or 3M Deutschland GmbH.

© 2016, 3M. All rights reserved.

ENGLISH

Product Description

Impregum™ Penta® Soft Quick is a fast-setting polyether impression material with a medium-bodied consistency. The mixing ratio, based on volume, is 5 parts base paste: 1 part catalyst. The paste is mixed in Pentamix versions from Pentamix™ 2 (September 1999) on, earlier Pentamix units (before September 1999) extrude the material too slow.

Use the mixing ratio to ensure the impression material is highly suitable for impressions with slight undercuts, even without additional blocking out. Each of the foil bags is sealed with a **PentaMatric™** cap. This **PentaMatric** cap opens the foil bag automatically in the Pentamix mixing unit as soon as enough pressure is built up by the plunger.

See the pertinent instructions for use for details on the Pentamix and accessories, Polyether Adhesive, Penta™ Elastomer Spritze and Penta™ Mixing Tips Red.

⚠ These instructions for use must be kept for reference for the duration of product use.

Indications

- Impressions of inlay, onlay, veneer, crown, and bridge preparations
- Fixation and implant impressions
- Impregum Penta Soft Quick is especially suitable for taking impressions of single-unit and two-unit preparations.

Precautionary Measures

For Patients

The base paste of this product contains substances that may cause an allergic reaction by skin contact in certain individuals. If prolonged contact with oral soft tissue occurs, flush with large amounts of water. If an adverse reaction is observed, a patient should discontinue use of the product in this patient and seek medical attention as needed. Known allergens are 1-Dodecylimidazole and mint flavors.

For Dental Personnel

The base paste of this product contains substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured product. If skin contact occurs, wash skin with soap and water. Use of protective gloves and a no-touch technique is recommended.

3M MSDSs can be obtained from www.mm.com or contact your local subsidiary.

Preparatory Steps

- Block-out undercuts or areas where gingival recession exists to prevent the impression material from “locking” onto the tooth structure. Failure to block-out may make tray removal difficult, or cause extraction of natural teeth or prosthesis. The block-out is especially important when using impression materials with high Shore hardness.

Impression tray:

- Impression trays generally used for precision impressions are suitable.
- For sufficient adhesion, apply a thin layer of Polyether Adhesive to the tray and allow it to dry completely according to the Polyether Adhesive Instructions for Use.

Pentamix/Penta Cartridge/Foil Bag

- Use the Impregum Penta Soft Quick foil bag only in combination with the designated Penta cartridge and with Penta Mixing Tips Red.
- Prepare new foil bags before taking the first impression by extruding about 3 cm of paste. The paste must emerge in a uniform color.
- If a new mixing tip has already been mounted when the cartridge is inserted, make sure that the drive shaft engages the mixing tip prior to start of mixing.

Retraction

Solutions based on aluminum hydroxide chloride or aluminum sulfate are suitable retraction agents for this product, with the exception of ephrine (adrenaline), 8-Hydroxyquinoline sulfate, and iron (III) sulfate. Iron (III) sulfate may impair the setting of polyether impressions materials.

- Keep the area dry where the impression is being taken.
- If necessary, use hemostatic threads for subgingival preparations.
- When taking the impression, remove completely any residue of the retraction agent by rinsing and drying.

Dosing and Mixing

- Dosing and mixing are done automatically in the Pentamix mixing unit.
- Only the Pentamix 2 or younger Pentamix versions must be used to mix Impregum Penta Soft Quick.

Times

Working Time from start of mixing	Initial setting Time	Setting Time from start of mixing
min:sec	min:sec	min:sec
Impregum Penta Soft Quick	Max. 1:00 + 3:00	= Max. 4:00

Symbols on product

Impregum Penta Soft Quick	Max. 1:00	+	3:00	=	Max. 4:00
---------------------------	-----------	---	------	---	-----------

The processing time (max. 1:00 min) begins with the entry of the paste into the Penta Mixing Tip Red (= start of mixing). It includes the filling of the Elastomer Spritze, application around the preparation, loading of the impression tray, and the positioning of the tray in the mouth.

The residence time in the mouth must always be 3 minutes, regardless of the previously noted processing time.

Impression Taking

Nono-phase technique

- To apply around the preparation, attach the Penta Elastomer Spritze to the Penta Mixing Tip Red of the Pentamix mixing unit and fill it.
- Then load the impression tray that has been prepared with adhesive. Keep the mixing tip immersed in the impression material all the times.
- Apply around the dried preparation from the bottom up. Keep the tip of the intra-oral tip immersed in the impression material and in contact with the tooth surface at all times during application.
- Application around the preparation and loading of the impression tray must be coordinated so that the loaded impression tray can be placed in the mouth immediately after the application.**
- The processing times shown for the materials used must not be exceeded.** Otherwise, the two materials will not set together. The consequences would be found in the impression.

- Position the loaded impression tray in the mouth immediately after the application and hold it in place, without applying pressure, until setting is complete.
- When the setting time has elapsed, remove the impression from the mouth. To remedy initial adhesion ("setting the valve"), especially with upper jaw impressions, remove the tray from the gingiva on one side in a posterior

position. If this proves difficult it may be necessary to carefully blow some air or water between the impression and the gingiva.

After Impression Taking

- Thoroughly examine and explore the sulcus of the prepared teeth and surrounding dentition. Remove any residual cured impression material from the mouth.

Hygiene

- Place the impression in a standard disinfectant solution commonly used for impressions for the period of time recommended by the manufacturer. Excessive disinfection may damage the impression.
- After disinfection, rinse the impression under running water for approx. 15 sec, then dry.

Model Preparation

- Prepare a cast from the impression with a commercial specialized stone plaster no earlier than 30 min and no later than 14 days after taking the impression.
- To avoid introducing bubbles into the model, briefly pre-rinse with water and dry with air. Do not use surfactants as these impair the quality of polyether impressions and are not necessary.
- Polyether impressions may be silver-coated, whereas copper-coating is not feasible.

Plating has not set may be removed with ethanol or by rinsing with water and soap.

- Acetone can be used to remove the adhesive from impression trays.

Notes

- All temperatures below 18° C/64° F: the viscosity of the pastes will increase to such an extent that there may be mixing problems in the unit. Store the pastes for 1 day at a minimum temperature of 18° C/64° F; they can then be processed without loss of quality.
- Direct exposure to sunlight and damp storage conditions may damage the impression, see Storage!
- Polyether impressions should not be exposed to solvent-containing liquids under any circumstances. Swelling and imprecise molding may result.
- 3M ESPE polyether materials may only be combined with 3M ESPE polyethers. The working and setting times of the individual products must be strictly observed.

Storage and Stability

- Store at 18–25° C/64–77° F. **Do not refrigerate!**
- Do not use after the expiration date.
- Store impressions in a dry place away from light at temperatures below 30° C/86° F.

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty

3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for its application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

Limitation of Liability

Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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